

Anti-Forkhead box protein P2 FoxP2 Monoclonal Antibody

Catalog Number: M01329


Overview

| | |
|----------------------|--|
| Product Name | Anti-Forkhead box protein P2 FoxP2 Monoclonal Antibody |
| Reactive Species | Human |
| Description | Boster Bio Anti-Forkhead box protein P2 FoxP2 Monoclonal Antibody catalog # M01329. Tested in ELISA, WB applications. This antibody reacts with Human. |
| Application | ELISA, WB |
| Clonality | Monoclonal 5C11A8 |
| Formulation | Ascitic fluid containing 0.03% sodium azide. |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles. |
| Host | Mouse |
| Uniprot ID | O15409 |

Technical Details

| | |
|---------------------|--|
| Immunogen | Purified recombinant fragment of human FoxP2 expressed in E. Coli. |
| Cross Reactivity | No cross reactivity with other proteins. |
| Isotype | IgG |
| Form | Liquid |
| Concentration | 1 mg/ml |
| Purification | Affinity purification |
| Suggested Dilutions | WB 1:500-1:2000 ELISA 1:10000 |

Anti-Forkhead box protein P2 FoxP2 Monoclonal Antibody (M01329) Images



Western Blot analysis using FoxP2 Monoclonal Antibody against recombinant protein.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Forkhead box protein P2 FoxP2 Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.