

Anti-GOLPH3 Antibody (C-term)

Catalog Number: M01333

About GOLPH3

Mediates the cis and medial Golgi localization of mannosyltransferases through direct binding of their cytosolic domains. Involved in modulation of mTOR signaling. Involved in the regulation of mitochondrial lipids, leading to increase of mitochondrial mass. Potential oncogene.

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview


Product Name	Anti-GOLPH3 Antibody (C-term)
Reactive Species	Human
Description	Boster Bio Anti-GOLPH3 Antibody (C-term) (Catalog # M01333). Tested in WB application(s). This antibody reacts with Human.
Application	WB
Clonality	Monoclonal 905CT9.1.1
Formulation	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Host	Mouse
Uniprot ID	Q9H4A6

Technical Details

Immunogen	Purified His-tagged GOLPH3 protein was used to produced this monoclonal antibody.
Predicted Reactive Species	Bovine, Chicken, Horse, Monkey, Pig, Sheep, Xenopus
Isotype	IgG1
Purification	This antibody is purified through a protein G column, followed by dialysis against PBS.
Suggested Dilutions	WB: 1:1000 Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions

	should be determined by end users.
--	------------------------------------

Anti-GOLPH3 Antibody (C-term) (M01333) Images



GOLPH3 Antibody (C-term) western blot analysis in 293,MDA-MB-231 cell line lysates (35ug/lane).This demonstrates the GOLPH3 antibody detected the GOLPH3 protein (arrow).

Anti-GOLPH3 Antibody (C-term)

For Research Use Only. Not for use in diagnostic procedures.