

Anti-c-Kit Rabbit Monoclonal Antibody

Catalog Number: M01335-1

About KIT

Activated by icilin, eucalyptol, menthol, cold and modulation of intracellular pH. Involved in menthol sensation. Permeable for monovalent cations sodium, potassium, and cesium and divalent cation calcium.

Overview


Product Name	Anti-c-Kit Rabbit Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-c-Kit Rabbit Monoclonal Antibody catalog # M01335-1. Tested in WB, IHC applications. This antibody reacts with Human.
Application	IHC, WB
Clonality	Monoclonal CGF-11
Formulation	Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol, 0.4-0.5mg/ml BSA.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P10721

Technical Details


Immunogen	A synthesized peptide derived from human c-Kit
Isotype	Rabbit IgG
Form	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Purification	Affinity-chromatography
Suggested Dilutions	Dilute the sample so that the expected range of concentrations fall within the detection range of this kit. If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples. Some PubMed article(s) citing the expression level of this target are as follows: Boster Bio's internal QC testing used: WB 1:500-1:2000

	IHC 1:50-1:200
--	----------------

Anti-c-Kit Rabbit Monoclonal Antibody (M01335-1) Images



Western blot analysis of c-Kit expression in Human fetal lung tissue lysate.



Immunohistochemical analysis of paraffin-embedded human breast cancer, using c-Kit Antibody .

1 Publications Citing This Product

1. PubMed ID: 28115023, Astaxanthin attenuates total body irradiation-induced hematopoietic system injury in mice via inhibition of oxidative stress and apoptosis

Visit bosterbio.com/anti-c-kit-rabbit-monoclonal-antibody-m01335-1-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-c-Kit Rabbit Monoclonal Antibody