

Anti-Hepcidin / HAMP Reference Antibody (Ludwig-Maximilians U. anti_Hepsin)

Catalog Number: M01347-3


Overview

Product Name	Anti-Hepcidin / HAMP Reference Antibody (Ludwig-Maximilians U. anti_Hepsin)
Reactive Species	Human
Description	Boster Bio Anti-Hepcidin / HAMP Reference Antibody (Ludwig-Maximilians U. anti_Hepsin) (Catalog # M01347-3). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Human. Endotoxin: < 0.570EU/ug,determined by LAL method. Expression system: CHO Cell
Application	ELISA, Flow Cytometry, Functional Assay, Kinetics
Clonality	Monoclonal
Formulation	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Storage Instructions	Maintain refrigerated at 2–8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.
Uniprot ID	P81172


Technical Details

Isotype	IgG1
Form	Liquid

Anti-Hepcidin / HAMP Reference Antibody (Ludwig-Maximilians U. anti_Hepsin) (M01347-3) Images



Anti-Hepcidin/HAMP Reference Antibody (Ludwig-Maximilians U. anti_Hepsin) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 90%



The purity of Anti-Hepcidin/HAMP Reference Antibody (Ludwig-Maximilians U. anti_Hepsin) is more than 95%

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Hepcidin / HAMP Reference Antibody (Ludwig-Maximilians U. anti_Hepsin)

For Research Use Only. Not for use in diagnostic procedures.