

## Anti-Cdc6 Rabbit Monoclonal Antibody

Catalog Number: M01355

### About CDC6

Furin is likely to represent the ubiquitous endoprotease activity within constitutive secretory pathways and capable of cleavage at the RX (K/R) R consensus motif.

### Overview

Product Name	Anti-Cdc6 Rabbit Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-Cdc6 Rabbit Monoclonal Antibody catalog # M01355. Tested in WB, IHC, ICC/IF applications. This antibody reacts with Human.
Conjugate	FITC
Application	IF, IHC, ICC, WB
Clonality	Monoclonal H1E-3
Formulation	Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol, 0.4-0.5mg/ml BSA.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q99741

### Technical Details

Immunogen	A synthesized peptide derived from human Cdc6
Predicted Reactive Species	Human, Primate
Cross Reactivity	Detects ~20kDa. Does not cross-react with alphaB-crystallin, betaL-crystallin, ̢H- crystallin, gamma-crystallin, HSP25, HSP27 or HSP47 proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Purification	Affinity-chromatography

**Suggested Dilutions**

Dilute the sample so that the expected range of concentrations fall within the detection range of this kit.

If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples.

Some PubMed article(s) citing the expression level of this target are as follows:


Boster Bio's internal QC testing used:

WB 1:500-1:2000

IHC 1:50-1:200

ICC/IF 1:50-1:200

## Anti-Cdc6 Rabbit Monoclonal Antibody (M01355) Images



Western blot analysis of Cdc6 expression in HeLa cell lysate treated with hydroxyurea.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Cdc6 Rabbit Monoclonal Antibody