

Anti-ER81 ETV1 Monoclonal Antibody

Catalog Number: M01383

Overview

Product Name	Anti-ER81 ETV1 Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-ER81 ETV1 Monoclonal Antibody catalog # M01383. Tested in ELISA, WB applications. This antibody reacts with Human.
Application	ELISA, WB
Clonality	Monoclonal 1C8B6
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	P50549

Technical Details

Immunogen	Purified recombinant fragment of ER81 (aa1-191) expressed in E. Coli.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Affinity purification
Suggested Dilutions	WB 1:500-1:2000 ELISA 1:10000

Anti-ER81 ETV1 Monoclonal Antibody (M01383) Images




Figure 1. Western blotting validation for Anti-ER81 ETV1 Monoclonal Antibody M01383

Western Blot (WB) analysis using ER81 Monoclonal Antibody against truncated Trx-ETV1 recombinant protein (1) and full-length ETV1 (aa1-477)-hlgGfc transfected CHO-K1 cell lysate(2).

Electrophoresis was performed on a SDS-PAGE gel. To determine SDS-PAGE gel concentration

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ER81 ETV1 Monoclonal Antibody