

Anti-ETV1 Antibody (C-Term)

Catalog Number: M01383-1

About ETV1

Transcriptional activator that binds to DNA sequences containing the consensus pentanucleotide 5'-CGGA[AT]-3'.


Overview

Product Name	Anti-ETV1 Antibody (C-Term)
Reactive Species	Human, Mouse
Description	Boster Bio Anti-ETV1 Antibody (C-Term) (Catalog # M01383-1). Tested in WB application(s). This antibody reacts with Human, Mouse.
Application	WB
Clonality	Polyclonal
Formulation	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P50549

Technical Details

Immunogen	This ETV1 antibody is generated from a rabbit immunized with a KLH conjugated synthetic peptide between 434-465 amino acids from human ETV1.
Predicted Reactive Species	Bovine, Chicken, Horse, Monkey, Pig, Sheep, Xenopus
Isotype	Rabbit IgG
Purification	This antibody is purified through a protein A column, followed by peptide affinity purification.
Suggested Dilutions	WB: 1:2000

Anti-ETV1 Antibody (C-Term) (M01383-1) Images



All lanes : Anti-ETV1 Antibody (C-Term) at 1:2000 dilution
Lane 1: HeLa whole cell lysates
Lane 2: HT-29 whole cell lysates
Lane 3: human testis lysates
Lane 4: mouse heart lysates
Lane 5: human brain lysates
Lane 6: human lung lysates
Lane 7: mouse kidney lysates
Lysates/proteins at 20 µg per lane. Secondary Goat Anti-Rabbit IgG, (H+L), Peroxidase conjugated at 1/10000 dilution
Predicted band size : 55 kDa
Blocking/Dilution buffer: 5% NFDm/TBST.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ETV1 Antibody (C-Term)

For Research Use Only. Not for use in diagnostic procedures.