

Anti-PKR EIF2AK2 Rabbit Monoclonal Antibody

Catalog Number: M01384

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Anti-PKR EIF2AK2 Rabbit Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-PKR EIF2AK2 Rabbit Monoclonal Antibody catalog # M01384. Tested in WB, IHC, ICC/IF, IP applications. This antibody reacts with Human.
Application	IP, IF, IHC, ICC, WB
Clonality	Monoclonal FCH-5
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P19525


Technical Details

Immunogen	A synthesized peptide derived from human PKR
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography


Suggested Dilutions

WB 1:1000-5000
IHC 1:50-200
ICC/IF 1:50-200
IP 1:20


Anti-PKR EIF2AK2 Rabbit Monoclonal Antibody (M01384) Images



Western blot analysis of PKR expression in MCF-7 cell lysate.



Immunohistochemical analysis of paraffin-embedded human kidney, using PKR Antibody.



Immunofluorescent analysis of HeLa cells, using PKR Antibody.

[Anti-PKR EIF2AK2 Rabbit Monoclonal Antibody](#)

For Research Use Only. Not for use in diagnostic procedures.