

Anti-Hexokinase-2 (N-term) monoclonal HK2 Antibody

Catalog Number: M01389


Overview

Product Name	Anti-Hexokinase-2 (N-term) monoclonal HK2 Antibody
Reactive Species	Human, Rat
Description	Boster Bio Anti-Hexokinase-2 (N-term) monoclonal HK2 Antibody catalog # M01389. Tested in WB,Flow Cytometry applications. This antibody reacts with Human,Rat.
Application	Flow Cytometry, WB
Clonality	Monoclonal 1C8B6
Formulation	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	P52789


Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids N-terminus of Human Hexokinase-2.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:1000 FC: 1:50-100

Anti-Hexokinase-2 (N-term) monoclonal HK2 Antibody (M01389) Images



Western blot detection of Hexokinase-2 (N-terminus) antibody in PC12, HeLa&Jurkat lysates using Hexokinase-2 (N-terminus) antibody (1:1000 diluted) .



Flow Cytometry analysis of K562 cells stained with Hexokinase-2(N-terminus) (red, 1/100 dilution), followed by FITC-conjugated goat anti-mouse IgG. Blue line histogram represents the isotype control, normal mouse IgG.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Hexokinase-2 (N-term) monoclonal HK2 Antibody

For Research Use Only. Not for use in diagnostic procedures.