

Anti-HXK II HK2 Monoclonal Antibody

Catalog Number: M01389-1


Overview

Product Name	Anti-HXK II HK2 Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-HXK II HK2 Monoclonal Antibody catalog # M01389-1. Tested in ELISA, Flow Cytometry, IHC, WB applications. This antibody reacts with Human.
Application	ELISA, Flow Cytometry, IF, IHC, WB
Clonality	Monoclonal 3D3
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	P52789


Technical Details

Immunogen	Purified recombinant fragment of human HXK II expressed in E. Coli.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Affinity purification
Suggested Dilutions	WB 1:500-1:2000 IHC 1:200-1:1000 FC 1:200-1:400 ELISA 1:10000 IF 1:50-200


Anti-HXK II HK2 Monoclonal Antibody (M01389-1) Images



Western Blot analysis using HXK II Monoclonal Antibody against Jurkat (1), HeLa (2) and HEK293 (3) cell lysate.



Immunohistochemistry analysis of paraffin-embedded esophagus cancer tissues (left) and human lung cancer (right) with DAB staining using HXK II Monoclonal Antibody.



Flow cytometric analysis of K562 cells using HXK II Monoclonal Antibody (green) and negative control (purple).

1 Publications Citing This Product

1. PubMed ID: 32768579, Chen J,Wu D,Dong Z,Chen A,Liu S.The expression and role of glycolysis-associated molecules in infantile hemangioma.Life Sci.2020 Oct 15;259:118215.doi:10.1016/j.lfs.2020.118215.Epub 2020 Aug 5.PMID:32768579.

Visit bosterbio.com/anti-hxk-ii-hk2-monoclonal-antibody-m01389-1-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-HXK II HK2 Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.