

## Anti-HXK II HK2 Monoclonal Antibody

Catalog Number: M01389-2


### Overview

Product Name	Anti-HXK II HK2 Monoclonal Antibody
Reactive Species	Human, Mouse, Pig, Rat
Description	Boster Bio Anti-HXK II HK2 Monoclonal Antibody catalog # M01389-2. Tested in WB applications. This antibody reacts with Human, Mouse, Rat, Swine.
Application	WB
Clonality	Monoclonal 3D3
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	P52789


### Technical Details

Immunogen	Purified recombinant human HXK II (N-terminus) protein fragments expressed in E.coli.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Affinity purification
Suggested Dilutions	WB 1:1000-1:2000

## Anti-HXK II HK2 Monoclonal Antibody (M01389-2) Images



Western Blot analysis using HXK II Monoclonal Antibody against PC12, HeLa, Jurkat cell lysate.



Flow cytometric analysis of K562 cells stained with HXK II Monoclonal Antibody (red), followed by FITC-conjugated goat anti-mouse IgG. Blue line histogram represents the isotype control, normal mouse IgG.

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-HXK II HK2 Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.