

Anti-IL-6Ra/CD126 Reference Antibody (tocilizumab)


Catalog Number: M01425-2

Overview


Product Name	Anti-IL-6Ra/CD126 Reference Antibody (tocilizumab)
Reactive Species	Cynomolgus monkey, Human
Description	Boster Bio Anti-IL-6Ra/CD126 Reference Antibody (tocilizumab) (Catalog # M01425-2). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Cynomolgus monkey, Human. Endotoxin: < 0.016EU/ug,determined by LAL method. Expression system: CHO Cell
Application	ELISA, Flow Cytometry, Functional Assay, Kinetics
Clonality	Monoclonal
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.
Uniprot ID	P08887

Technical Details


Anti-IL-6Ra/CD126 Reference Antibody (tocilizumab) (M01425-2) Images




Anti-IL-6Ra/CD126 Reference Antibody (tocilizumab) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95%



The purity of Anti-IL-6Ra/CD126 Reference Antibody (tocilizumab) is more than 95%



Human IL-6R CHO cells were stained with Anti-IL-6Ra/CD126 Reference Antibody (tocilizumab) and negative control protein respectively



Anti-IL-6Ra/CD126 Reference Antibody (tocilizumab)-induced FACS Blocking activity was evaluated using Human IL-6R CHO cells. The IC50 was approximately 1.303nM.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-IL-6Ra/CD126 Reference Antibody (tocilizumab)

For Research Use Only. Not for use in diagnostic procedures.