

Anti-Arrestin-beta-2 ARRB2 Monoclonal Antibody

Catalog Number: M01443


Overview

Product Name	Anti-Arrestin-beta-2 ARRB2 Monoclonal Antibody
Reactive Species	Human, Mouse, Pig, Rat
Description	Boster Bio Anti-Arrestin-beta-2 ARRB2 Monoclonal Antibody catalog # M01443. Tested in IF, WB applications. This antibody reacts with Human, Mouse, Rat, Swine.
Application	IF, WB
Clonality	Monoclonal 3H6D12 / 4H9C8 / 4H9F1
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	P32121


Technical Details

Immunogen	Purified recombinant human Arrestin-beta-2 (C-terminus) protein fragments expressed in E.coli.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Affinity purification
Suggested Dilutions	WB 1:1000-1:2000 IF 1:100-1:500

Anti-Arrestin-beta-2 ARRB2 Monoclonal Antibody (M01443) Images



Western Blot analysis using Arrestin-beta-2 Monoclonal Antibody against HepG2 cell lysate.



Immunofluorescence analysis of HeLa cells using Arrestin-beta-2 Monoclonal Antibody.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Arrestin-beta-2 ARRB2 Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.