

## Anti-ACVRL1 / ALK-1 Reference Antibody (ascrinvacumab)

Catalog Number: M01468-2


### Overview

Product Name	Anti-ACVRL1 / ALK-1 Reference Antibody (ascrinvacumab)
Reactive Species	Cynomolgus monkey, Human
Description	Boster Bio Anti-ACVRL1 / ALK-1 Reference Antibody (ascrinvacumab) (Catalog # M01468-2). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Cynomolgus monkey, Human. Endotoxin: < 0.586EU/ug,determined by LAL method. Expression system: CHO Cell
Application	ELISA, Flow Cytometry, Functional Assay, Kinetics
Clonality	Monoclonal
Formulation	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.
Uniprot ID	P37023


### Technical Details

Isotype	IgG2SA
Form	Liquid

## Anti-ACVRL1 / ALK-1 Reference Antibody (ascrivacumab) (M01468-2) Images



Anti-ACVRL1/ALK-1 Reference Antibody (ascrivacumab) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 90%



The purity of Anti-ACVRL1/ALK-1 Reference Antibody (ascrivacumab) is more than 95%

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ACVRL1 / ALK-1 Reference Antibody (ascrivacumab)

For Research Use Only. Not for use in diagnostic procedures.