

## Anti-Ephrin type-A receptor 5 EphA5 Monoclonal Antibody

Catalog Number: M01494


### Overview

Product Name	Anti-Ephrin type-A receptor 5 EphA5 Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-Ephrin type-A receptor 5 EphA5 Monoclonal Antibody catalog # M01494. Tested in ELISA, WB applications. This antibody reacts with Human.
Application	ELISA, WB
Clonality	Monoclonal 8B10B1 / 8B10F5
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	P54756

### Technical Details

Immunogen	Purified recombinant fragment of EphA5 (aa620-774) expressed in E. Coli.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Affinity purification
Suggested Dilutions	WB 1:500-1:2000 ELISA 1:10000

## Anti-Ephrin type-A receptor 5 EphA5 Monoclonal Antibody (M01494) Images



Western blotting validation for Anti-Ephrin type-A receptor 5 EphA5 Monoclonal Antibody M01494 Western Blot (WB) analysis using EphA5 Monoclonal Antibody against truncated EPHA5-His recombinant protein (1) and truncated EPHA5(aa620-774)-hIgGFc transfected CHO-K1 cell lysate(2). Electrophoresis was performed on a SDS-PAGE gel. To determine SDS-PAGE gel concentration

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Ephrin type-A receptor 5 EphA5 Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.