

Anti-HMBS Mouse Monoclonal Antibody [Clone ID: OTI3F4]

Catalog Number: M01506-1


Overview

Product Name	Anti-HMBS Mouse Monoclonal Antibody [Clone ID: OTI3F4]
Reactive Species	Human, Mouse, Rat
Description	Boster Bio HMBS mouse monoclonal antibody, clone OTI3F4 (formerly 3F4). Catalog# M01506-1. Tested in IHC, WB. This antibody reacts with Human, Mouse, Rat.
Application	IHC, WB
Clonality	Monoclonal OTI3F4
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	P08397


Technical Details

Immunogen	Full length human recombinant protein of human HMBS (NP_000181) produced in HEK293T cell.
Isotype	IgG1
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	WB 1:2000


Anti-HMBS Mouse Monoclonal Antibody [Clone ID: OTI3F4] (M01506-1) Images




HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY HMBS (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-HMBS.




Immunohistochemical staining of paraffin-embedded Adenocarcinoma of Human ovary tissue using anti-HMBS mouse monoclonal antibody. (M01506-1) Dilution: 1:150




Immunohistochemical staining of paraffin-embedded Carcinoma of Human pancreas tissue using anti-HMBS mouse monoclonal antibody. (M01506-1) Dilution: 1:150



Immunohistochemical staining of paraffin-embedded Human lymphoma tissue using anti-HMBS mouse monoclonal antibody. (M01506-1) Dilution: 1:150



Immunohistochemical staining of paraffin-embedded Human Kidney tissue within the normal limits using anti-HMBS mouse monoclonal antibody. (M01506-1) Dilution: 1:150



Immunohistochemical staining of paraffin-embedded Carcinoma of Human kidney tissue using anti-HMBS mouse monoclonal antibody. (M01506-1) Dilution: 1:150

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-HMBS Mouse Monoclonal Antibody [Clone ID: OTI3F4]

For Research Use Only. Not for use in diagnostic procedures.