

## Anti-Siglec-3 / CD33 Reference Antibody (Gemtuzumab ozogamicin)

Catalog Number: M01508-6


### Overview

Product Name	Anti-Siglec-3 / CD33 Reference Antibody (Gemtuzumab ozogamicin)
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Siglec-3 / CD33 Reference Antibody (Gemtuzumab ozogamicin) (Catalog # M01508-6). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Rat, Human, Mouse. Endotoxin: < 0.599EU/ug,determined by LAL method. Expression system: CHO Cell
Application	ELISA, Flow Cytometry, Functional Assay, Kinetics
Clonality	Monoclonal
Formulation	Purified monoclonal antibody supplied in 100mM Pro-Ac, 20mM Arg, pH5.0, without preservative.This antibody is purified through a protein A column.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.
Uniprot ID	P20138


### Technical Details

Isotype	IgG1
Form	Liquid

## Anti-Siglec-3 / CD33 Reference Antibody (Gemtuzumab ozogamicin) (M01508-6) Images



Anti-Siglec-3/CD33 Reference Antibody (Gemtuzumab ozogamicin) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 90%



The purity of Anti-Siglec-3/CD33 Reference Antibody (Gemtuzumab ozogamicin) is more than 95%

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Siglec-3 / CD33 Reference Antibody (Gemtuzumab ozogamicin)

For Research Use Only. Not for use in diagnostic procedures.