

Anti-ROBO1 Reference Antibody (Asclepius Technology patent anti-Robo1 CAR)

Catalog Number: M01530

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.




Overview

Product Name	Anti-ROBO1 Reference Antibody (Asclepius Technology patent anti-Robo1 CAR)
Reactive Species	Human
Description	Boster Bio Anti-ROBO1 Reference Antibody (Asclepius Technology patent anti-Robo1 CAR) (Catalog # M01530). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Human. Endotoxin: < 0.605EU/ug,determined by LAL method. Expression system: CHO Cell
Application	ELISA, Flow Cytometry, Functional Assay, Kinetics
Clonality	Monoclonal
Formulation	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Storage Instructions	Maintain refrigerated at 2–8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.
Uniprot ID	Q9Y6N7


Technical Details

Isotype	IgG1
Form	Liquid

Anti-ROBO1 Reference Antibody (Asclepius Technology patent anti-Robo1 CAR) (M01530) Images



Anti-ROBO1 Reference Antibody (Asclepius Technology patent anti-Robo1 CAR) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95%



The purity of Anti-ROBO1 Reference Antibody (Asclepius Technology patent anti-Robo1 CAR) is more than 95%

Anti-ROBO1 Reference Antibody (Asclepius Technology patent anti-Robo1 CAR)

For Research Use Only. Not for use in diagnostic procedures.