

Anti-MSPR / RON / CD136 Reference Antibody (Aveo anti-RON)

Catalog Number: M01532-4

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.




Overview

Product Name	Anti-MSPR / RON / CD136 Reference Antibody (Aveo anti-RON)
Reactive Species	Human
Description	Boster Bio Anti-MSPR / RON / CD136 Reference Antibody (Aveo anti-RON) (Catalog # M01532-4). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Human. Endotoxin: < 0.608EU/ug,determined by LAL method. Expression system: CHO Cell
Application	ELISA, Flow Cytometry, Functional Assay, Kinetics
Clonality	Monoclonal
Formulation	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.
Uniprot ID	Q04912


Technical Details

Isotype	IgG1
Form	Liquid


Anti-MSPR / RON / CD136 Reference Antibody (Aveo anti-RON) (M01532-4) Images



Anti-MSPR/RON/CD136 Reference Antibody (Aveo anti-RON) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95%



The purity of Anti-MSPR/RON/CD136 Reference Antibody (Aveo anti-RON) is more than 97.3%



Hs746T cells were stained with Anti-MSPR/RON/CD136 Reference Antibody (Aveo anti-RON) and negative control protein respectively

[Anti-MSPR / RON / CD136 Reference Antibody \(Aveo anti-RON\)](#)

For Research Use Only. Not for use in diagnostic procedures.