

Anti-Bag3 Rabbit Monoclonal Antibody

Catalog Number: M01547-2

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.




Overview

Product Name	Anti-Bag3 Rabbit Monoclonal Antibody
Reactive Species	Human, Mouse, Rat
Description	This Recombinant Rabbit Monoclonal Antibody detects endogenous levels of Bag3 protein. Validated for research with WB,IHC,IF,IP,ELISA applications. Anti Bag3 is reactive for Human, Mouse, Rat samples
Application	ELISA, IP, IF, IHC, WB
Clonality	Monoclonal
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05% stabilizing protein This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C/1 year
Host	Rabbit
Uniprot ID	O95817


Technical Details

Isotype	IgG, Kappa
Purification	Protein A
Suggested Dilutions	IHC 1:200-1:1000 WB 1:1000-1:5000 IF 1:200-1:1000 ELISA 1:5000-1:20000 IP 1:50-1:200


Anti-Bag3 Rabbit Monoclonal Antibody (M01547-2) Images




Human breast carcinoma was stained with anti-Bag3 rabbit antibody



Mouse lung was stained with anti-Bag3 rabbit antibody



Rat lung was stained with anti-Bag3 rabbit antibody



Various whole cell lysates were separated by 4-20% SDS-PAGE, and the membrane was blotted with anti-Bag3 antibody. The HRP-conjugated Goat anti-Rabbit IgG(H + L) antibody was used to detect the antibody. Lane 1: HEK293 Lane 2: HeLa Lane 3: C6 Lane 4: 3T3-L1 Predicted band size: 62kDa Observed band size: 80kDa

Anti-Bag3 Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.