

## Anti-Integrin alpha V ITGAV Rabbit Monoclonal Antibody

Catalog Number: M01561


### Overview

Product Name	Anti-Integrin alpha V ITGAV Rabbit Monoclonal Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Integrin alpha V ITGAV Rabbit Monoclonal Antibody catalog # M01561. Tested in WB, IHC, ICC/IF, IP applications. This antibody reacts with Human, Mouse, Rat.
Application	IP, IF, IHC, ICC, WB
Clonality	Monoclonal ABGC-9
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P06756


### Technical Details

Immunogen	A synthesized peptide derived from human Integrin alpha V
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography
Suggested Dilutions	WB 1:500-2000 IHC 1:50-200 ICC/IF 1:50-200 IP 1:20


## Anti-Integrin alpha V ITGAV Rabbit Monoclonal Antibody (M01561) Images



Western blot analysis of ITGAV using anti-ITGAV antibody (M01561). Electrophoresis was performed on a 5-20% SDS-PAGE gel at 70V (Stacking gel) / 90V (Resolving gel) for 2-3 hours. The sample well of each lane was loaded with 30 ug of sample under reducing conditions. Lane 1: human A549 whole cell lysates, Lane 2: human A431 whole cell lysates, Lane 3: human 293T whole cell lysates, Lane 4: human MCF-7 whole cell lysates, Lane 5: rat spleen tissue lysates, Lane 6: rat brain tissue lysates, Lane 7: mouse spleen tissue lysates, Lane 8: mouse brain tissue lysates. After electrophoresis, proteins were transferred to a nitrocellulose membrane at 150 mA for 50-90 minutes. Blocked the membrane with 5% non-fat milk/TBS for 1.5 hour at RT. The membrane was incubated with rabbit anti-ITGAV antigen affinity purified monoclonal antibody (Catalog # M01561) at 1:500 overnight at 4°C, then washed with TBS-0.1%Tween 3 times with 5 minutes each and probed with a goat anti-rabbit IgG-HRP secondary antibody at a dilution of 1:1000 for 1.5 hour at RT. The signal is developed using an Enhanced Chemiluminescent detection (ECL) kit (Catalog # EK1002) with Tanon 5200 system. A specific band was detected for ITGAV at approximately 130, 140, 150 kDa. The expected band size for ITGAV is at 116 kDa.





Immunofluorescent analysis using the Antibody at 1:50 dilution.



Immunofluorescent analysis using the Antibody at 1:50 dilution.

Immunofluorescent analysis using the Antibody at 1:150 dilution.



Immunofluorescent analysis using the Antibody at 1:500 dilution.

## 4 Publications Citing This Product

1. PubMed ID: 25666351, Chen Q, Zhang Z, Liu J, He Q, Zhou Y, Shao G, Sun X, Cao X, Gong A, Jiang P. Mol Cells. 2015 Mar;38(3):221-8. Doi: 10.14348/Molcells.2015.2170. Epub 2015 Feb 4. A Fibrin Matrix Promotes The Differentiation Of Emscs Isolated From Nasal Respiratory ...
2. PubMed ID: 29312560, Xuan J, Chen Y, Zhu L, Guo Y, Deng L, Zheng Y, Wang Z, Wang Z, Ao M. Oncotarget. 2017 Sep 28;8(65):108676-108691. doi: 10.18632/oncotarget.21358. eCollection 2017 Dec 12. Ultrasound molecular imaging with cRGD-PLGA-PFOB nanoparticles for liver fib...
3. PubMed ID: 21747684, Study on the Expression and Clinical Significances of Lewis y Antigen and Integrin  $\alpha$ v,  $\beta$ 3 in Epithelial Ovarian Tumors

Visit [bosterbio.com/anti-integrin-alpha-v-rabbit-monoclonal-antibody-m01561-boster.html](http://bosterbio.com/anti-integrin-alpha-v-rabbit-monoclonal-antibody-m01561-boster.html) to see all 4 publications.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Integrin alpha V ITGAV Rabbit Monoclonal Antibody  
For Research Use Only. Not for use in diagnostic procedures.