

Anti-Integrin α V β 3 (ITGAV & ITGB3) Reference Antibody (etaracizumab)

Catalog Number: M01561-3


Overview

Product Name	Anti-Integrin α V β 3 (ITGAV & ITGB3) Reference Antibody (etaracizumab)
Reactive Species	Human
Description	Boster Bio Anti-Integrin α V β 3 (ITGAV & ITGB3) Reference Antibody (etaracizumab) (Catalog # M01561-3). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Human. Endotoxin: < 0.610EU/ug,determined by LAL method. Expression system: CHO Cell
Application	ELISA, Flow Cytometry, Functional Assay, Kinetics
Clonality	Monoclonal
Formulation	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Storage Instructions	Maintain refrigerated at 2–8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.
Uniprot ID	P06756 & P05106


Technical Details

Isotype	IgG1
Form	Liquid

Anti-Integrin α V β 3 (ITGAV & ITGB3) Reference Antibody (etaracizumab) (M01561-3) Images



Anti-Integrin α V β 3 (ITGAV & ITGB3) Reference Antibody (etaracizumab) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 90%



The purity of Anti-Integrin α V β 3 (ITGAV & ITGB3) Reference Antibody (etaracizumab) is more than 95%

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Integrin α V β 3 (ITGAV & ITGB3) Reference Antibody (etaracizumab)

For Research Use Only. Not for use in diagnostic procedures.