

Anti-BubR1/BUB1B Antibody Picoband® (monoclonal, 8B3) HRP Conjugated

Catalog Number: M01564-3-HRP

About BUB1B

Mitotic checkpoint serine/threonine-protein kinase BUB1 beta is an enzyme that in humans is encoded by the BUB1B gene. This gene encodes a kinase involved in spindle checkpoint function. The protein has been localized to the kinetochore and plays a role in the inhibition of the anaphase-promoting complex/cyclosome (APC/C), delaying the onset of anaphase and ensuring proper chromosome segregation. Impaired spindle checkpoint function has been found in many forms of cancer.

Overview

Product Name	Anti-BubR1/BUB1B Antibody Picoband® (monoclonal, 8B3) HRP Conjugated
Reactive Species	Human, Rat
Clonality	Monoclonal 8B3
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Mouse
Uniprot ID	O60566

Technical Details

Immunogen	E.coli-derived human BubR1/BUB1B recombinant protein (Position: K26-E448).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Mouse IgG2b
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-BubR1/BUB1B Antibody (monoclonal, 8B3) - HRP

For Research Use Only. Not for use in diagnostic procedures.