

## Anti-AIF/AIFM1 Antibody Picoband® (monoclonal, 2I5) Fluoro488 Conjugated

Catalog Number: M01571-1-Fluoro488

### About AIFM1

Apoptosis-inducing factor 1, mitochondrial, also known as AIF or PDCD8 is a protein that in humans is encoded by the AIFM1 gene. AIFM1 gene is mapped to Xq26.1 based on an alignment of the AIFM1 sequence with the genomic sequence. This gene encodes a flavoprotein essential for nuclear disassembly in apoptotic cells, and it is found in the mitochondrial intermembrane space in healthy cells. Induction of apoptosis results in the translocation of this protein to the nucleus where it affects chromosome condensation and fragmentation. In addition, this gene product induces mitochondria to release the apoptogenic proteins cytochrome c and caspase-9. Mutations in this gene cause combined oxidative phosphorylation deficiency 6, which results in a severe mitochondrial encephalomyopathy. A related pseudogene has been identified on chromosome 10.

### Overview

Product Name	Anti-AIF/AIFM1 Antibody Picoband® (monoclonal, 2I5) Fluoro488 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (Flow Cytometry, IF, IHC, ICC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Monoclonal 2I5
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% Na <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Mouse
Uniprot ID	O95831

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human AIF, identical to the related mouse and rat sequences.
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Mouse IgG1
Form	Liquid
Concentration	0.5 mg/mL
Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm

Suggested Dilutions

Optimal dilutions should be determined by end users.

### 3 Publications Citing This Product

1. PubMed ID: 27895670, Wnt5a Increases Properties of Lung Cancer Stem Cells and Resistance to Cisplatin through Activation of Wnt5a/PKC Signaling Pathway
2. PubMed ID: 26645545, Toll-like receptor 4 contributes to chronic itch, allodynia, and spinal astrocyte activation in male mice
3. PubMed ID: 27330750, Mechanism of apoptosis induction in human hepatocellular carcinoma cells following treatment with a gecko peptides mixture

Visit [bosterbio.com/anti-aif-picoband-trade-antibody-monoclonal-m01571-1-boster.html](http://bosterbio.com/anti-aif-picoband-trade-antibody-monoclonal-m01571-1-boster.html) to see all 3 publications.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-AIF/AIFM1 Antibody (monoclonal, 215) - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.