

Anti-N-Cadherin Rabbit Monoclonal Antibody

Catalog Number: M01577-6

Overview

Product Name	Anti-N-Cadherin Rabbit Monoclonal Antibody
Reactive Species	Human, Mouse
Description	This Recombinant Rabbit Monoclonal Antibody detects endogenous levels of N-Cadherin (ABT156R) rabbit mAb protein. Validated for research with WB, IHC, ELISA applications. Anti N-Cadherin is reactive for Human, Mouse, (predicted: Rat) samples
Application	ELISA, IHC, WB
Clonality	Monoclonal
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05% stabilizing protein This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C/1 year
Host	Rabbit
Uniprot ID	P19022

Technical Details

Predicted Reactive Species	Rat
Isotype	IgG, Kappa
Purification	Protein A

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-N-Cadherin Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.