

Anti-CDA Antibody (Center)

Catalog Number: M01578

About CDA

This enzyme scavenges exogenous and endogenous cytidine and 2'-deoxycytidine for UMP synthesis.


Overview

Product Name	Anti-CDA Antibody (Center)
Reactive Species	Human
Description	Boster Bio Anti-CDA Antibody (Center) (Catalog # M01578). Tested in WB application(s). This antibody reacts with Human.
Application	WB
Clonality	Monoclonal 921CT14.2.1
Formulation	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Host	Mouse
Uniprot ID	P32320

Technical Details

Immunogen	Purified His-tagged CDA protein was used to produce this monoclonal antibody.
Predicted Reactive Species	Bovine, Chicken, Horse, Monkey, Pig, Sheep, Xenopus
Isotype	IgG1
Purification	This antibody is purified through a protein G column, followed by dialysis against PBS.
Suggested Dilutions	WB: 1:1000

Anti-CDA Antibody (Center) (M01578) Images



CDA Antibody (Center) western blot analysis in HeLa cell line lysates (35ug/lane). This demonstrates the CDA antibody detected the CDA protein (arrow).

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CDA Antibody (Center)

For Research Use Only. Not for use in diagnostic procedures.