

Anti-GDF15 Antibody

Catalog Number: M01583

About GDF15

Regulates food intake, energy expenditure and body weight in response to metabolic and toxin-induced stresses (PubMed:[28953886](#), PubMed:[28846097](#), PubMed:[28846098](#), PubMed:[28846099](#), PubMed:[23468844](#), PubMed:[29046435](#)). Binds to its receptor, GFRAL, and activates GFRAL-expressing neurons localized in the area postrema and nucleus tractus solitarius of the brainstem (PubMed:[28953886](#), PubMed:[28846097](#), PubMed:[28846098](#), PubMed:[28846099](#)). It then triggers the activation of neurons localized within the parabrachial nucleus and central amygdala, which constitutes part of the 'emergency circuit' that shapes feeding responses to stressful conditions (PubMed:[28953886](#)). On hepatocytes, inhibits growth hormone signaling (By similarity).


Overview

Product Name	Anti-GDF15 Antibody
Reactive Species	Human
Description	Boster Bio Anti-GDF15 Antibody (Catalog # M01583). Tested in WB application(s). This antibody reacts with Human.
Application	WB
Clonality	Monoclonal 2182CT228.18.8
Formulation	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Host	Mouse
Uniprot ID	Q99988

Technical Details

Immunogen	This GDF15 antibody is generated from a mouse immunized with a recombinant protein from the human GDF15.
Predicted Reactive Species	Human
Isotype	IgG1,kappa
Purification	This antibody is purified through a protein G column, followed by dialysis against PBS.
Suggested Dilutions	WB: 1:4000

Anti-GDF15 Antibody (M01583) Images



Anti-GDF15 Antibody at 1:4000 dilution + HT-1080 whole cell lysate Lysates/proteins at 20 μ g per lane. Secondary Goat Anti-mouse IgG, (H+L), Peroxidase conjugated at 1/10000 dilution. Predicted band size : 34 kDa Blocking/Dilution buffer: 5% NFDM/TBST.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-GDF15 Antibody

For Research Use Only. Not for use in diagnostic procedures.