

Anti-(Mouse) Nr4a2 Antibody (Center)

Catalog Number: M01592

About Nr4a2

This enzyme scavenge exogenous and endogenous cytidine and 2'-deoxycytidine for UMP synthesis.

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.




Overview

Product Name	Anti-(Mouse) Nr4a2 Antibody (Center)
Reactive Species	Mouse, Rat
Description	Boster Bio Anti-(Mouse) Nr4a2 Antibody (Center) (Catalog # M01592). Tested in WB, IHC-P application(s). This antibody reacts with Mouse, Rat.
Application	IHC-P, WB
Clonality	Polyclonal
Formulation	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q06219


Technical Details

Immunogen	This mouse Nr4a2 antibody is generated from a rabbit immunized with a KLH conjugated synthetic peptide between 164-197 amino acids from the Central region of mouse Nr4a2.
Predicted Reactive Species	Bovine, Chicken, Horse, Monkey, Pig, Sheep, Xenopus
Isotype	Rabbit IgG
Purification	This antibody is purified through a protein A column, followed by peptide affinity purification.
Suggested Dilutions	WB: 1:2000 IHC-P: 1:25


Anti-(Mouse) Nr4a2 Antibody (Center) (M01592) Images




All lanes : Anti-Nr4a2 Antibody (Center) at 1:2000 dilution
Lane 1: mouse brain lysates
Lane 2: rat brain lysates
Lysates/proteins at 20 µg per lane. Secondary Goat Anti-Rabbit IgG, (H+L), Peroxidase conjugated at 1/10000 dilution
Predicted band size : 67 kDa
Blocking/Dilution buffer: 5% NFDm/TBST.



Anti-Nr4a2 Antibody (Center) at 1:2000 dilution + mouse pancreas lysates
Lysates/proteins at 20 µg per lane. Secondary Goat Anti-Rabbit IgG, (H+L), Peroxidase conjugated at 1/10000 dilution
Predicted band size : 67 kDa
Blocking/Dilution buffer: 5% NFDm/TBST.



M01592 staining (Mouse) Nr4a2 in mouse brain sections by Immunohistochemistry (IHC-P -paraformaldehyde-fixed, paraffin-embedded sections). Tissue was fixed with formaldehyde and blocked with 3% BSA for 0.5 hour at room temperature; antigen retrieval was by heat mediation with a citrate buffer (pH6). Samples were incubated with primary antibody (1/25) for 1 hours at 37°C. A undiluted biotinylated goat polyvalent antibody was used as the secondary antibody.



M01592 staining (Mouse) Nr4a2 in mouse spleen sections by Immunohistochemistry (IHC-P -paraformaldehyde-fixed, paraffin-embedded sections). Tissue was fixed with formaldehyde and blocked with 3% BSA for 0.5 hour at room temperature; antigen retrieval was by heat mediation with a citrate buffer (pH6). Samples were incubated with primary antibody (1/25) for 1 hours at 37°C. A undiluted biotinylated goat polyvalent antibody was used as the secondary antibody.

[Anti-\(Mouse\) Nr4a2 Antibody \(Center\)](#)

For Research Use Only. Not for use in diagnostic procedures.