

Anti-eNOS NOS3 Monoclonal Antibody

Catalog Number: M01604


Overview

Product Name	Anti-eNOS NOS3 Monoclonal Antibody
Reactive Species	Human, Mouse, Rabbit, Rat
Description	Boster Bio Anti-eNOS NOS3 Monoclonal Antibody catalog # M01604. Tested in WB applications. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Monoclonal Mix
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	P29474

Technical Details

Immunogen	Recombinant Protein of eNOS
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Suggested Dilutions	WB 1:500-2000

Anti-eNOS NOS3 Monoclonal Antibody (M01604) Images



Western blot analysis of lysates from 1) Rat Heart Tissue, 2) huvec, 3) Jurkat cells, (Green primary antibody was diluted at 1:1000, 4° over night, secondary antibody was diluted at 1:10000, 37° 1hour. (Red Actin beta Polyclonal Antibody was diluted at 1:5000 as loading control, 4° over night, secondary antibody was diluted at 1:10000, 37° 1hour.

1 Publications Citing This Product

1. PubMed ID: 31309562, Cai P, Lu Z, Jiang T, Wang Z, Yang Y, Zheng L, Zhao J. Syndecan-4 involves in the pathogenesis of rheumatoid arthritis by regulating the inflammatory response and apoptosis of fibroblast-like synoviocytes. *J Cell Physiol*. 2020 Feb;235(2):1746-1758. doi:10.1002/jcp

Visit bosterbio.com/anti-enos-nos3-monoclonal-antibody-m01604-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-eNOS NOS3 Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.