

## Anti-Cytokeratin 1 Rabbit Monoclonal Antibody

Catalog Number: M01639


### Overview

Product Name	Anti-Cytokeratin 1 Rabbit Monoclonal Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Cytokeratin 1 Rabbit Monoclonal Antibody catalog # M01639. Tested in WB application. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Monoclonal GHO-11
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P04264

### Technical Details

Immunogen	A synthesized peptide derived from human Cytokeratin 1
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography
Suggested Dilutions	WB 1:5000-1:20000

## Anti-Cytokeratin 1 Rabbit Monoclonal Antibody (M01639) Images



Western blot analysis of Cytokeratin 1 expression in Human skin lysate.

## 2 Publications Citing This Product

1. PubMed ID: 26870216, A metanephric adenoma of the kidney associated with polycythemia: A case report
2. PubMed ID: 24962416, Stanniocalcin-1 protects bovine intestinal epithelial cells from oxidative stress-induced damage

Visit [bosterbio.com/anti-cytokeratin-1-rabbit-monoclonal-antibody-m01639-boster.html](https://bosterbio.com/anti-cytokeratin-1-rabbit-monoclonal-antibody-m01639-boster.html) to see all 2 publications.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Cytokeratin 1 Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.