

## Anti-MUC16/Ca125 Rabbit Monoclonal Antibody

Catalog Number: M01641


### Overview

Product Name	Anti-MUC16/Ca125 Rabbit Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-MUC16/Ca125 Rabbit Monoclonal Antibody catalog # M01641. Tested in WB, ICC/IF, Flow Cytometry applications. This antibody reacts with Human.
Application	Flow Cytometry, IF, ICC, WB
Clonality	Monoclonal GOG-13
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q8WXI7

### Technical Details

Immunogen	A synthesized peptide derived from human MUC16
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography
Suggested Dilutions	WB 1:1000-5000 ICC/IF 1:50-200 FC 1:20

## Anti-MUC16/Ca125 Rabbit Monoclonal Antibody (M01641) Images



Western blot analysis of MUC16 expression in Human ovary cancer lysate.

### 1 Publications Citing This Product

1. PubMed ID: 32765681, Wang H, Yan C, Ye H. Overexpression of MUC16 predicts favourable prognosis in MUC16-mutant cervical cancer related to immune response. *Exp Ther Med*. 2020 Aug;20(2):1725-1733. doi:10.3892/etm.2020.8836. Epub 2020 Jun 4. PMID:32765681; PMCID:PMC7388571.

Visit [bosterbio.com/anti-muc16-rabbit-monoclonal-antibody-m01641-boster.html](https://bosterbio.com/anti-muc16-rabbit-monoclonal-antibody-m01641-boster.html) to see all 1 publications.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-MUC16/Ca125 Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.