

Anti-ALP Reference Antibody (Seagen Patent Anti-Alpp)

Catalog Number: M01718-5


Overview

Product Name	Anti-ALP Reference Antibody (Seagen Patent Anti-Alpp)
Reactive Species	Cynomolgus monkey, Human
Description	Boster Bio Anti-ALP Reference Antibody (Seagen Patent Anti-Alpp) (Catalog # M01718-5). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Cynomolgus monkey, Human. Endotoxin: < 0.633EU/ug,determined by LAL method. Expression system: CHO Cell
Application	ELISA, Flow Cytometry, Functional Assay, Kinetics
Clonality	Monoclonal
Formulation	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.
Uniprot ID	P05187-1


Technical Details

Isotype	IgG1
Form	Liquid


Anti-ALP Reference Antibody (Seagen Patent Anti-Alpp) (M01718-5) Images



Anti-ALP Reference Antibody (Seagen Patent Anti-Alpp) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95%



The purity of Anti-ALP Reference Antibody (Seagen Patent Anti-Alpp) is more than 97.51%



Immobilized human ALPP His at 4 μg/mL can bind Anti-ALP Reference Antibody (Seagen Patent Anti-Alpp)

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ALP Reference Antibody (Seagen Patent Anti-Alpp)

For Research Use Only. Not for use in diagnostic procedures.