

Anti-DOCK8 Mouse Monoclonal Antibody [Clone ID: OTI8D9]

Catalog Number: M01771-2


Overview

Product Name	Anti-DOCK8 Mouse Monoclonal Antibody [Clone ID: OTI8D9]
Reactive Species	Human, Mouse, Rat
Description	Boster Bio DOCK8 mouse monoclonal antibody, clone OTI8D9 (formerly 8D9). Catalog# M01771-2. Tested in FC, WB. This antibody reacts with Human, Mouse, Rat.
Application	Flow Cytometry, WB
Clonality	Monoclonal OTI8D9
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	Q8NF50


Technical Details

Immunogen	Human recombinant protein fragment corresponding to amino acids 84-405 of human DOCK8 (NP_982272) produced in E.coli.
Isotype	IgG1
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	WB 1:500

Anti-DOCK8 Mouse Monoclonal Antibody [Clone ID: OTI8D9] (M01771-2) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY DOCK8 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-DOCK8 (1:500).



Flow cytometric Analysis of penetrated k562 cells

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-DOCK8 Mouse Monoclonal Antibody [Clone ID: OTI8D9]

For Research Use Only. Not for use in diagnostic procedures.