

## Anti-CD166 ALCAM Monoclonal Antibody

Catalog Number: M01788-1


### Overview

Product Name	Anti-CD166 ALCAM Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-CD166 ALCAM Monoclonal Antibody catalog # M01788-1. Tested in ELISA, IHC applications. This antibody reacts with Human.
Application	ELISA, IF, IHC
Clonality	Monoclonal 10F1G12
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	Q13740

### Technical Details

Immunogen	Purified recombinant fragment of CD166 (aa405-524) expressed in E. Coli.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Affinity purification
Suggested Dilutions	IHC 1:200-1:1000 ELISA 1:10000 IF 1:50-200

## Anti-CD166 ALCAM Monoclonal Antibody (M01788-1) Images



Immunohistochemistry validation of ALCAM using Anti-CD166 ALCAM Monoclonal Antibody (M01788-1). Immunohistochemistry (IHC) analysis of paraffin-embedded Human Ovary carcinoma (A)

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CD166 ALCAM Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.