

## Anti-Acetyl Coenzyme A Carboxylase Rabbit Monoclonal Antibody

Catalog Number: M01802-2

### Overview

Product Name	Anti-Acetyl Coenzyme A Carboxylase Rabbit Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-Acetyl Coenzyme A Carboxylase Rabbit Monoclonal Antibody catalog # M01802-2. Tested in WB, IHC, ICC/IF applications. This antibody reacts with Human.
Application	IF, IHC, ICC, WB
Clonality	Monoclonal 32A06
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q13085

### Technical Details

Immunogen	A synthesized peptide derived from human Acetyl Coenzyme A Carboxylase
Isotype	IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography
Suggested Dilutions	WB 1:500-2000 IHC 1:50-200 ICC/IF 1:50-200

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Acetyl Coenzyme A Carboxylase Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.