

Anti-Thromboxane A2 receptor Rabbit Monoclonal Antibody

Catalog Number: M01812


Overview

Product Name	Anti-Thromboxane A2 receptor Rabbit Monoclonal Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Thromboxane A2 receptor Rabbit Monoclonal Antibody catalog # M01812. Tested in WB application. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Monoclonal 27T87
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P21731

Technical Details

Immunogen	A synthesized peptide derived from human Thromboxane A2 receptor
Isotype	IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography
Suggested Dilutions	WB 1:500-2000

Anti-Thromboxane A2 receptor Rabbit Monoclonal Antibody (M01812) Images



Western blot analysis of Thromboxane A2 receptor expression in human platelet cell lysate.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Thromboxane A2 receptor Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.