

Anti-Spermine synthase Rabbit Monoclonal Antibody

Catalog Number: M01831-1

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview


Product Name	Anti-Spermine synthase Rabbit Monoclonal Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Spermine synthase Rabbit Monoclonal Antibody catalog # M01831-1. Tested in WB, ICC/IF applications. This antibody reacts with Human, Mouse, Rat.
Application	IF, ICC, WB
Clonality	Monoclonal 27S72
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P52788

Technical Details


Immunogen	A synthesized peptide derived from human Spermine synthase
Isotype	IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography

Suggested Dilutions	WB 1:500-2000 ICC/IF 1:50-200
---------------------	----------------------------------


Anti-Spermine synthase Rabbit Monoclonal Antibody (M01831-1) Images




All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



Western blot analysis of Spermine synthase expression in HeLa cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

[Anti-Spermine synthase Rabbit Monoclonal Antibody](#)

For Research Use Only. Not for use in diagnostic procedures.