

## Anti-PER3 Rabbit Monoclonal Antibody

Catalog Number: M01835


### Overview

Product Name	Anti-PER3 Rabbit Monoclonal Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-PER3 Rabbit Monoclonal Antibody catalog # M01835. Tested in WB applications. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Monoclonal 24P84
Formulation	Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol, 0.4-0.5mg/ml BSA.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P56645

### Technical Details

Immunogen	A synthesized peptide derived from human PER3
Isotype	IgG
Form	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Purification	Affinity-chromatography
Suggested Dilutions	Dilute the sample so that the expected range of concentrations fall within the detection range of this kit. If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples. Some PubMed article(s) citing the expression level of this target are as follows: Boster Bio's internal QC testing used WB 1:500-1:2000  For protocols, please visit <a href="https://www.bosterbio.com/protocol-and-troubleshooting/">https://www.bosterbio.com/protocol-and-troubleshooting/</a>

### Anti-PER3 Rabbit Monoclonal Antibody (M01835) Images



Western blot analysis of PER3 expression in K562 cell lysate.

## 1 Publications Citing This Product

1. PubMed ID: 10.1080/13880209.2021.1961823, Kaempferol-induced GPER upregulation attenuates atherosclerosis via the PI3K/AKT/Nrf2 pathway

Visit [bosterbio.com/anti-per3-rabbit-monoclonal-antibody-m01835-boster.html](https://bosterbio.com/anti-per3-rabbit-monoclonal-antibody-m01835-boster.html) to see all 1 publications.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PER3 Rabbit Monoclonal Antibody