

## Anti-Apolipoprotein L1 ApoL1 Monoclonal Antibody

Catalog Number: M01841


### Overview

Product Name	Anti-Apolipoprotein L1 ApoL1 Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-Apolipoprotein L1 ApoL1 Monoclonal Antibody catalog # M01841. Tested in ELISA, WB applications. This antibody reacts with Human.
Application	ELISA, WB
Clonality	Monoclonal 1D4
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	O14791

### Technical Details

Immunogen	Purified recombinant fragment of ApoL1 expressed in E. Coli.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Affinity purification
Suggested Dilutions	WB 1:500-1:2000 ELISA 1:10000

## Anti-Apolipoprotein L1 ApoL1 Monoclonal Antibody (M01841) Images



Western blotting validation for Anti-Apolipoprotein L1 ApoL1 Monoclonal Antibody M01841 Western Blot (WB) analysis using ApoL1 Monoclonal Antibody against human plasma (1). Electrophoresis was performed on a SDS-PAGE gel. To determine SDS-PAGE gel concentration

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Apolipoprotein L1 ApoL1 Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.