

Anti-PRDX1/Peroxiredoxin 1 Rabbit Monoclonal Antibody

Catalog Number: M01845

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Anti-PRDX1/Peroxiredoxin 1 Rabbit Monoclonal Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-PRDX1/Peroxiredoxin 1 Rabbit Monoclonal Antibody catalog # M01845. Tested in WB, IHC, ICC/IF, Flow Cytometry applications. This antibody reacts with Human, Mouse.
Application	Flow Cytometry, IF, IHC, ICC, WB
Clonality	Monoclonal IIG-16
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q06830


Technical Details

Immunogen	A synthesized peptide derived from human PRDX1
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography

Suggested Dilutions

WB 1:500-2000
IHC 1:50-200
ICC/IF 1:50-200
FC 1:20

Anti-PRDX1/Peroxiredoxin 1 Rabbit Monoclonal Antibody (M01845) Images



Western blot analysis of PRDX1/Peroxiredoxin 1 using anti-PRDX1/Peroxiredoxin 1 antibody (M01845). Electrophoresis was performed on a 5-20% SDS-PAGE gel at 70V (Stacking gel) / 90V (Resolving gel) for 2-3 hours. The sample well of each lane was loaded with 30 ug of sample under reducing conditions. Lane 1: human A549 whole cell lysates, Lane 2: human Hela whole cell lysates, Lane 3: human HepG2 whole cell lysates, Lane 4: human U251 whole cell lysates, Lane 5: mouse liver tissue lysates, Lane 6: mouse spleen tissue lysates. After electrophoresis, proteins were transferred to a nitrocellulose membrane at 150 mA for 50-90 minutes. Blocked the membrane with 5% non-fat milk/TBS for 1.5 hour at RT. The membrane was incubated with rabbit anti-PRDX1/Peroxiredoxin 1 antigen affinity purified monoclonal antibody (Catalog # M01845) at 1:1000 overnight at 4°C, then washed with TBS-0.1%Tween 3 times with 5 minutes each and probed with a goat anti-rabbit IgG-HRP secondary antibody at a dilution of 1:500 for 1.5 hour at RT. The signal is developed using an Enhanced Chemiluminescent detection (ECL) kit (Catalog # EK1002) with Tanon 5200 system. A specific band was detected for PRDX1/Peroxiredoxin 1 at approximately 24 kDa. The expected band size for PRDX1/Peroxiredoxin 1 is at 22,24 kDa.

Anti-PRDX1/Peroxiredoxin 1 Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.