

Anti-SCDGFB (PDGFD) Mouse Monoclonal Antibody [Clone ID: OTI1C1]

Catalog Number: M01875


Overview

Product Name	Anti-SCDGFB (PDGFD) Mouse Monoclonal Antibody [Clone ID: OTI1C1]
Reactive Species	Human, Mouse, Rat
Description	Boster Bio PDGFD mouse monoclonal antibody, clone OTI1C1 (formerly 1C1). Catalog# M01875. Tested in FC, WB. This antibody reacts with Human, Mouse, Rat.
Application	Flow Cytometry, WB
Clonality	Monoclonal OTI1C1
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	Q9GZP0


Technical Details

Immunogen	Human recombinant protein fragment corresponding to amino acids 19-247 of human PDGFD (NP_079484) produced in E.coli.
Isotype	IgG1
Concentration	0.97 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	WB: 1:2000 Flow cytometry: 1:100

Anti-SCDGFB (PDGFD) Mouse Monoclonal Antibody [Clone ID: OTI1C1] (M01875) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY PDGFD (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-PDGFD.



HEK293T cells transfected with either PDGFD (Myc-DDK-tagged) overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-PDGFD antibody (M01875)

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-SCDGFB (PDGFD) Mouse Monoclonal Antibody [Clone ID: OTI1C1]

For Research Use Only. Not for use in diagnostic procedures.