

Anti-c-Rel Monoclonal Antibody

Catalog Number: M01880-2

Overview

Product Name	Anti-c-Rel Monoclonal Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-c-Rel Monoclonal Antibody catalog # M01880-2. Tested in ELISA, IF, IHC, WB applications. This antibody reacts with Human, Mouse.
Application	ELISA, IF, WB
Clonality	Monoclonal 1E7
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	Q04864

Technical Details

Immunogen	Purified recombinant fragment of human c-Rel expressed in E. Coli.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Affinity purification
Suggested Dilutions	WB 1:500-1:2000 IF 1:200-1:1000 ELISA 1:10000

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-c-Rel Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.