

Anti-Alpha A Crystallin Antibody Picoband® (monoclonal, 10B9) HRP Conjugated

Catalog Number: M01900-2-HRP

About CRYAA

Alpha-crystallin A chain is a protein that in humans is encoded by the CRYAA gene. Mammalian lens crystallins are divided into alpha, beta, and gamma families. Alpha crystallins are composed of two gene products: alpha-A and alpha-B, for acidic and basic, respectively. Alpha crystallins can be induced by heat shock and are members of the small heat shock protein (HSP20) family. They act as molecular chaperones although they do not renature proteins and release them in the fashion of a true chaperone; instead they hold them in large soluble aggregates. Two additional functions of alpha crystallins are an autokinase activity and participation in the intracellular architecture. The encoded protein has been identified as a moonlighting protein based on its ability to perform mechanistically distinct functions. Alpha-A and alpha-B gene products are differentially expressed; alpha-A is preferentially restricted to the lens and alpha-B is expressed widely in many tissues and organs. Defects in this gene cause autosomal dominant congenital cataract (ADCC).

Overview

| | |
|----------------------|---|
| Product Name | Anti-Alpha A Crystallin Antibody Picoband® (monoclonal, 10B9) HRP Conjugated |
| Reactive Species | Human, Mouse, Rat |
| Clonality | Monoclonal 10B9 |
| Formulation | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ . |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. |
| Host | Mouse |
| Uniprot ID | P02489 |

Technical Details

| | |
|------------------|---|
| Immunogen | E. coli-derived human Alpha A Crystallin recombinant protein (Position: M1-S173). Human Alpha A Crystallin shares 94.8% amino acid (aa) sequence identity with both mouse and rat Alpha A Crystallin. |
| Cross Reactivity | No cross-reactivity with other proteins. |
| Isotype | Mouse IgG1 |
| Form | Liquid |
| Concentration | 0.5 mg/mL |
| Purification | Immunogen affinity purified. |
| Conjugate | HRP |

Suggested Dilutions

The intended application should be selected according to the customer's experimental requirements.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Alpha A Crystallin Antibody (monoclonal, 10B9) - HRP

For Research Use Only. Not for use in diagnostic procedures.