

Anti-Integrin alpha 5 ITGA5 Rabbit Monoclonal Antibody

Catalog Number: M01911


Overview

Product Name	Anti-Integrin alpha 5 ITGA5 Rabbit Monoclonal Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Integrin alpha 5 ITGA5 Rabbit Monoclonal Antibody catalog # M01911. Tested in WB, IHC, ICC/IF, IP, Flow Cytometry applications. This antibody reacts with Human, Mouse, Rat.
Application	Flow Cytometry, IP, IF, IHC, ICC, WB
Clonality	Monoclonal IDB-9
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P08648


Technical Details

Immunogen	A synthesized peptide derived from human Integrin alpha 5
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography
Suggested Dilutions	WB 1:500-2000 IHC 1:50-200 ICC/IF 1:50-200 IP 1:20 FC 1:20

Anti-Integrin alpha 5 ITGA5 Rabbit Monoclonal Antibody (M01911) Images



Western blot analysis of Integrin alpha 5 expression in Mouse heart lysate.



Immunofluorescent analysis of MCF-7 cells, using Integrin alpha 5 Antibody.

3 Publications Citing This Product

1. PubMed ID: 19918321, Characterization of Rat Hair Follicle Stem Cells Selected by Vario Magnetic Activated Cell Sorting System
2. PubMed ID: 26680220, Silencing of Paternally Expressed Gene 10 Inhibits Trophoblast Proliferation and Invasion
3. PubMed ID: 24643246, A Small Physiological Electric Field Mediated Responses of Extravillous Trophoblasts Derived from HTR8/SVneo Cells: Involvement of Activation of Focal Adhesion Kinase Signaling

Visit bosterbio.com/anti-integrin-alpha-5-rabbit-monoclonal-antibody-m01911-boster.html to see all 3 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Integrin alpha 5 ITGA5 Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.