

Anti-GPC3 / Glypican-3 Reference Antibody (Codrituzumab-MMAE)

Catalog Number: M01922-3

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.




Overview

Product Name	Anti-GPC3 / Glypican-3 Reference Antibody (Codrituzumab-MMAE)
Reactive Species	Human
Description	Boster Bio Anti-GPC3 / Glypican-3 Reference Antibody (Codrituzumab-MMAE) (Catalog # M01922-3). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Human. Endotoxin: < 0.660EU/ug,determined by LAL method. Expression system: CHO Cell
Application	ELISA, Flow Cytometry, Functional Assay, Kinetics
Clonality	Monoclonal
Formulation	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.
Uniprot ID	P51654


Technical Details

Isotype	IgG1
Form	Liquid


Anti-GPC3 / Glypican-3 Reference Antibody (Codrituzumab-MMAE) (M01922-3) Images



Anti-GPC3/Glypican-3 Reference Antibody (Codrituzumab-MMAE) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 90%



The purity of Anti-GPC3/Glypican-3 Reference Antibody (Codrituzumab-MMAE) is more than 99.73%



The killing ratio Codrituzumab-MMAE by HuH-7 increased with the increase of antibody concentration

Anti-GPC3 / Glypican-3 Reference Antibody (Codrituzumab-MMAE)

For Research Use Only. Not for use in diagnostic procedures.