

## Anti-RASA1 Mouse Monoclonal Antibody [Clone ID: OTI1H2]

Catalog Number: M01929


### Overview

Product Name	Anti-RASA1 Mouse Monoclonal Antibody [Clone ID: OTI1H2]
Reactive Species	Human, Mouse, Rat
Description	Boster Bio RASA1 (GAP ) mouse monoclonal antibody, clone OTI1H2 (formerly 1H2). Catalog# M01929. Tested in IF, WB. This antibody reacts with Human, Mouse, Rat.
Application	IF, WB
Clonality	Monoclonal OTI1H2
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	P20936


### Technical Details

Immunogen	Full length human recombinant rotein of human RASA1 (NP_002881) produced in HEK293T cell.
Isotype	IgG1
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	WB: 1:2000 IF: 1:100

## Anti-RASA1 Mouse Monoclonal Antibody [Clone ID: OTI1H2] (M01929) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY RASA1 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-RASA1.



Anti-RASA1 mouse monoclonal antibody (M01929) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY RASA1.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-RASA1 Mouse Monoclonal Antibody [Clone ID: OTI1H2]

For Research Use Only. Not for use in diagnostic procedures.