

## Anti-Integrin a2b3 (ITGA2 & ITGB3) Reference Antibody (tadocizumab)

Catalog Number: M01933-1


### Overview

Product Name	Anti-Integrin a2b3 (ITGA2 & ITGB3) Reference Antibody (tadocizumab)
Reactive Species	Human
Description	Boster Bio Anti-Integrin a2b3 (ITGA2 & ITGB3) Reference Antibody (tadocizumab) (Catalog # M01933-1). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Human. Endotoxin: < 0.661EU/ug,determined by LAL method. Expression system: CHO Cell
Application	ELISA, Flow Cytometry, Functional Assay, Kinetics
Clonality	Monoclonal
Formulation	Purified monoclonal antibody supplied in 0.1M Glycine,0.2M NaCl,(tris),pH6.0, without preservative.This antibody is purified through a protein A column.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.
Uniprot ID	P17301 & P05106


### Technical Details

Isotype	IgG1
Form	Liquid

## Anti-Integrin a2b3 (ITGA2 & ITGB3) Reference Antibody (tadocizumab) (M01933-1) Images



Anti-Integrin a2b3 (ITGA2 & ITGB3) Reference Antibody (tadocizumab) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 90%



The purity of Anti-Integrin a2b3 (ITGA2 & ITGB3) Reference Antibody (tadocizumab) is more than 95%

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Integrin a2b3 (ITGA2 & ITGB3) Reference Antibody (tadocizumab)

For Research Use Only. Not for use in diagnostic procedures.