

Anti-Calpain 1 Antibody Picoband® (monoclonal, 2I10) Fluoro550 Conjugated

Catalog Number: M01943-3-Fluoro550

About CAPN1

CAPN1 is also known as CANP or muCL. The calpains, calcium-activated neutral proteases, are nonlysosomal, intracellular cysteine proteases. The mammalian calpains include ubiquitous, stomach-specific, and muscle-specific proteins. The ubiquitous enzymes consist of heterodimers with distinct large, catalytic subunits associated with a common small, regulatory subunit. This gene encodes the large subunit of the ubiquitous enzyme, calpain 1. Several transcript variants encoding two different isoforms have been found for this gene.

Overview

| | |
|----------------------|--|
| Product Name | Anti-Calpain 1 Antibody Picoband® (monoclonal, 2I10) Fluoro550 Conjugated |
| Reactive Species | Human |
| Application | Flow Cytometry |
| Clonality | Monoclonal 2I10 |
| Formulation | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ . |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light. |
| Host | Mouse |
| Uniprot ID | P07384 |

Technical Details

| | |
|---------------------|---|
| Immunogen | E.coli-derived human Calpain 1 recombinant protein (Position: Q396-A555). Human Calpain 1 shares 86% amino acid (aa) sequence identity with both mouse and rat Calpain 1. |
| Cross Reactivity | No cross-reactivity with other proteins. |
| Isotype | Mouse IgG2a |
| Form | Liquid |
| Concentration | 0.5 mg/mL |
| Purification | Immunogen affinity purified. |
| Conjugate | Fluoro550 Excitation Wavelength: 562 nm Emission Wavelength: 576 nm |
| Suggested Dilutions | Flow Cytometry, Optimal dilutions should be determined by end users. |

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Calpain 1 Antibody (monoclonal, 2110) - Fluoro550

For Research Use Only. Not for use in diagnostic procedures.