

Anti-Desmin Antibody Picoband® (monoclonal, 2B5) Biotin Conjugated

Catalog Number: M01948-3-Biotin

Overview

| | |
|----------------------|--|
| Product Name | Anti-Desmin Antibody Picoband® (monoclonal, 2B5) Biotin Conjugated |
| Reactive Species | Human, Mouse, Rat |
| Application | WB, IHC, ELISA |
| Clonality | Monoclonal 2B5 |
| Formulation | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% Na ₃ N. |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. |
| Host | Mouse |
| Uniprot ID | P17661 |

Technical Details

| | |
|---------------------|---|
| Immunogen | E.coli-derived human Desmin recombinant protein (Position: M1-T304). Human Desmin shares 97% amino acid (aa) sequence identity with both mouse and rat Desmin. |
| Cross Reactivity | No cross-reactivity with other proteins. |
| Isotype | Mouse IgG2b |
| Form | Liquid |
| Concentration | 0.5 mg/mL |
| Conjugate | Biotin |
| Suggested Dilutions | Western blot, Optimal dilutions should be determined by end users. Immunohistochemistry (Paraffin-embedded Section), Optimal dilutions should be determined by end users. ELISA, Optimal dilutions should be determined by end users. |

7 Publications Citing This Product

1. PubMed ID: 10.1186/s12967-015-0603-0, Identification of apoptosis-related microRNAs and their target genes in myocardial infarction post-transplantation with skeletal myoblasts
2. PubMed ID: 10.1016/j.brainres.2015.08.003, Intranasal administration of human umbilical cord mesenchymal stem cells-conditioned medium

enhances vascular remodeling after stroke

3. PubMed ID: 10.1016/j.cellbi.2008.03.023, 17beta-Estradiol affects proliferation and apoptosis of rat prostatic smooth muscle cells by modulating cell cycle transition and related proteins

Visit bosterbio.com/anti-desmin-picoband-trade-antibody-monoclonal-m01948-3-boster.html to see all 7 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Desmin Antibody (monoclonal, 2B5) - Biotin

For Research Use Only. Not for use in diagnostic procedures.