

## Anti-Desmin Antibody Picoband® (monoclonal, 2B5) Fluoro647 Conjugated

Catalog Number: M01948-3-Fluoro647

### Overview

Product Name	Anti-Desmin Antibody Picoband® (monoclonal, 2B5) Fluoro647 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Flow Cytometry
Clonality	Monoclonal 2B5
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Mouse
Uniprot ID	P17661

### Technical Details

Immunogen	E.coli-derived human Desmin recombinant protein (Position: M1-T304). Human Desmin shares 97% amino acid (aa) sequence identity with both mouse and rat Desmin.
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Mouse IgG2b
Form	Liquid
Concentration	0.5 mg/mL
Conjugate	Fluoro647 Excitation Wavelength: 650 nm Emission Wavelength: 665 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

### 7 Publications Citing This Product

1. PubMed ID: 10.1186/s12967-015-0603-0, Identification of apoptosis-related microRNAs and their target genes in myocardial infarction post-transplantation with skeletal myoblasts
2. PubMed ID: 10.1016/j.brainres.2015.08.003, Intranasal administration of human umbilical cord mesenchymal stem cells-conditioned medium

enhances vascular remodeling after stroke

3. PubMed ID: 10.1016/j.cellbi.2008.03.023, 17beta-Estradiol affects proliferation and apoptosis of rat prostatic smooth muscle cells by modulating cell cycle transition and related proteins

Visit [bosterbio.com/anti-desmin-picoband-trade-antibody-monoclonal-m01948-3-boster.html](http://bosterbio.com/anti-desmin-picoband-trade-antibody-monoclonal-m01948-3-boster.html) to see all 7 publications.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Desmin Antibody (monoclonal, 2B5) - Fluoro647

For Research Use Only. Not for use in diagnostic procedures.